Fore more information, please visit our website at:
www.itouchrelief.com

MD-9

User Manual

EN

ES

User Manual

English

Contents

Foreword	2	Unique benefits	15
General information	3	Regular application principles	16
Overview	4	Regular application methods	17
Charging the MD-9 battery	5	Pain in the trunk	17
Features and mode description	6	Pain in the joints and limbs	18
Directions for use	7	Electrode placement variations	19
Electrode guidelines	8	Troubleshooting	20
Electrode recommendations	9	Contraindications	22
Maintenance and care	9	Warnings	23
Environmental specifications	10	Precautions	25
Technical specifications	10	Adverse reactions	25
Pulse parameters	11	Location of stimulation	25
Mode description, TENS	12	Conformity to safety standards	26
Mode description, EMS	13	Warranty information	27
Indications of use	14		

Foreword

Before operation, please read this user manual carefully and ensure you understand all instructions.

Do you suffer from **PAIN, MUSCLE TENSION**, or **STRESS**?

The ITouch Relief MD-9 can quickly relax tense muscles and relieve pain. It is a modern home therapy device that uses the latest bio-electrical technology.

The ITouch Relief MD-9 is designed to deliver an endless amount of relaxing, deep, soothing sensations exactly where you need relief. Relax, feel good, and rejuvenate your body and mind. Don't continue living with unhealthy physical and mental tension.

Enjoy the MD-9 and live life to the fullest; whether at home, at work, or on vacation. You'll love how good it feels!

Please read the instructions carefully before using this device.

The MD-9 is a dual-channel battery-operated muscle stimulation system designed to help relieve minor muscular aches and pains. It consists of two main components: an electronic stimulation module that generates the necessary signals and two sets of pro-neo neonatal ECG electrode pads with lead wires.

It is intended for temporary pain relief associated with sore and aching muscles in the lower back, upper extremities (arms), lower extremities (legs), and joints caused by exercise or routine household activities.

It is also intended to stimulate healthy muscles to improve and support muscle performance.

General Information

Included in this package:

MD-9 unit
1 set of LARGE single sided adhesive electrode pads (2.4 inches - round shape)
1 set of EXTRA-LARGE single sided adhesive electrode pads (3¾ inch - rectangular shape)
2 Electrode Wires
USB Cable
AC Adapter
User Manual
Pad Holder


Structure:

- | | |
|--------------------|--------------------------|
| 1. Control unit | 7. Output |
| 2. USB cable | 8. LCD screen |
| 3. AC adapter | 9. Power switch |
| 4. Electrode wires | 10. USB port |
| 5. Electrode pads | 11. A/B Channel selector |
| 6. Pad Holder | |

* In case of any discrepancy between the physical product and the images shown, please refer to the actual product.

Overview: control unit and accesories

USB Charger: Used to charge the MD-9 battery.


AC Adapter: Used with the USB cable to charge the battery.

Electrodes: Attach to bare skin to perform stimulation.

Electrode wires: Connect the device to the electrodes.

Pad holder: Used to store the electrode pads and wires after use.

Charging the MD-9 Battery

* Via a Standard Wall Outlet

- 1 Plug the AC adapter into any standard wall outlet.


- 2 Connect the small end of the USB cable to the MD-9 unit and the larger end to the adapter.

- 3 Charging is complete when the battery icon on the screen indicates full.

* Via a Computer USB Port

- 1 Connect the small end of the USB cable to the MD-9 unit and the larger end to a USB port on a computer (ensure the computer is powered on).

- 2 Charging is complete when the battery icon on the screen indicates full.

* Notes:

1. The lithium battery can be recharged using either the AC adapter or the USB cable.
2. **While charging:** If the unit is ON, the screen will display the charging progress. If the unit is OFF, the screen will remain blank. However, the battery will charge in either the ON or OFF position.
3. Charge the battery for 8–10 hours before first use.

Features

FDA approved Class II Medical Device

- Large, easy-to-read LCD display
- Backlight for increased readability
- Adjustable timer (10–60 minutes)
- 2-year warranty
- 8 stimulation modes
- Dual output channels
- Rechargeable internal lithium battery
- 20 intensity level settings

Mode Description

We recommend experimenting with each of the 8 modes to determine which provides the most comfortable and effective sensation for your current condition. You may find that a particular mode or combination of modes works best for you.

- Set the intensity to a level where you feel gentle muscle vibration and mild involuntary movements. The sensation should be comfortable, never painful.
- Each session should last at least 10 to 20 minutes.
- Use the device 1 to 3 times per day, depending on your needs.

Directions for Use

Three Key Points for Effective Operation:

1. **Target the Pain Point:** Apply the electrode pads directly on the area where pain is most commonly felt.
2. **Intensity:** The intensity must be gradually increased and it's better to reach the highest level you can comfortably tolerate.
3. **Application Duration:** Use the device 2–3 times per day, for 40–60 minutes per session, over a 10-day course.

- Connect the electrode wire to the output port on the top of the MD-9 unit.
- Snap the electrode pads onto the wire.
- Remove the protective film and place both pads on the skin, ensuring they are on muscle areas and not overlapping.
- Turn on the unit using the On/Off switch on top of the device.
- When you have turned on the unit, it will automatically start in Mode 1.
- To switch modes, press the M button.
- Once you have selected a mode, gradually increase the intensity using the “+” button. To reduce intensity, press the “–” button.
- To set the session time, press the T button.
- Sit back and enjoy the deep, soothing sensations!

* Notes:

The MD-9 unit is designed with your safety in mind. The output intensity can only be increased by pressing the “+” button, ensuring gradual and controlled adjustment. Even at the maximum intensity level, the output remains within a safe range. When switching modes, the device will automatically reset the intensity to the minimum level, allowing you to readjust comfortably for the new mode.

Electrode guidelines

- 1 Use only the electrodes provided with your MD-9 device. Third-party electrodes may have incompatible electrical properties and could pose a safety risk.
- 2 Do not share electrodes between users. Each person should use their own set to avoid skin reactions or the risk of cross-contamination.
- 3 Always turn off the device before removing or repositioning the electrodes.
- 4 Thoroughly wash and dry the skin before applying the electrodes to ensure proper adhesion and performance.
- 5 Apply the entire surface of each electrode firmly and evenly to the skin. Do not use electrodes that do not stick properly or only partially adhere.
- 6 If you notice skin redness under the electrodes after a session, wait until the skin has fully recovered before reapplying them to the same area.
- 7 Using electrode pads that are too small or improperly positioned may cause discomfort or even skin burns.

Electrode Recommendations

- The adhesive electrodes are supplied as a standard accessory with your MD-9 unit. These pads are long-lasting and reusable.
- The recommended electrode pad sizes are bigger round shape pads for use on large areas such as the shoulder, waist, abdomen, and leg muscles; and small oval pads for more compact areas like the forearm muscles, joints, etc.

Maintenance and Care

- 1 Ensure your skin is clean and free from dirt, oil, or lotions.
- 2 Before applying the pads, it is recommended to spray ITouch Relief conductor on your fingers and gently rub it on both pads. This helps maintain the pads' adhesiveness.
- 3 After using the device, turn it off.
- 4 Remove the pads carefully.
- 5 Reapply the protective films to the pads or place the pads on the pad holder.
- 6 Disconnect the wire from the device and return it to the storage box.
- 7 Store the MD-9 unit in a cool, dry place until next use.

- 9 -

Environmental Specifications

Operating conditions:

Temperature: 5°C to 40°C (41°F to 104°F)

Humidity (non-condensing): 30% to 75%

Atmospheric pressure: 700 to 1,060 HPa

Storage and transportation conditions:

Temperature: -10°C to 55°C (14°F to 131°F)

Humidity (non-condensing): 10% to 90%

Atmospheric pressure: 700 to 1,060 HPa

Technical Specifications

Dimensions: 86mm x 52mm x 11mm

Net weight: 52g

Power supply: DC 3.7V

- 10 -

Pulse Parameters

MD-9	TABLE: PULSE parameters			
Rated Supply Voltage (V)		DC 3.7 V		
Load Resistance (Ω)		1000 Ω		
PULSE DURATION(μ s)	PULSE repetition frequencies (Hz)	PULSE repetition amplitudes (V)	DC component (V)	Remarks
100	54	50~90	--	Mode 1
100	6	50~90	--	Mode 2
100	1.28	50~90	--	Mode 3
100	54	50~90	--	Mode 4
100	64	50~90	--	Mode 5
100	55	50~90	--	Mode 6
100	54	50~90	--	Mode 7.1; 7.2
100	<71	50~90	--	Mode 8.1; 8.2; 8.3; 8.4;
100	<60	50~90	--	Mode 9.1; 9.2; 9.3;
100	<48	50~90	--	Mode 10.1; 10.2;
100	8.7	50~90	--	Mode 11
100	<63	50~90	--	Mode 12.1; 12.2; 12.3; 12.4;

- 11 -

Mode description when used as TENS (Pain Relief)

Mode 3: Tapping – Simulates a secondary heartbeat. Tapping the muscle increases the body's natural blood flow and oxygen supply to the target area. It also helps flush out inflammation or swelling in the targeted area.

Mode 4: Hot Stone – Relaxation mode. Simulates a manual rolling massage to help the soft tissue release and relax.

Mode 6: Shuffle – A combination mode, ideal for cooling down after exercise.

Mode 8: Ramping Wave – Begins with Mode 3 (Tapping) and gradually increases speed and intensity, penetrating deeper into the muscle until it reaches a full contraction. Then it releases the muscle. This mode supports endurance and fatigue resistance.

Mode 10: Reflexology – Designed for use with the shoes. Simulates a reflexology session through socks or shoes using varied stimulation patterns.

Mode 11: Auricular Therapy – Also known as "alpha stim" or "microcurrent." Use with the earflap attachment. It simulates alpha brain waves and may support sleep patterns, behavior regulation, and endorphin release. It also relaxes the suboccipital muscles and promotes vasodilation, helping relieve migraines and headaches.

- 12 -

Mode description when used as EMS (Muscle Stimulation)

Mode 1: Kneading – Simulates the thumbs, fingers, and palms of a real massage therapist. The sensation typically feels like a vibration. This mode targets hypertension and stress, helping to relieve tension and restore full range of motion and flexibility.

Mode 2: Acupressure – Delivers a rapid, focused pulse aimed at muscle knots or adhesions. It applies enough pressure to break them down and release any built-up acidic (lactic) content within the muscle tissue.

Mode 5: Cupping – Simulates a grabbing and stretching motion on the muscle. It helps to lengthen, broaden, and then release the muscle, improving flexibility and mobility.

Mode 7: Contract and Hold – The longest contraction cycle on the device. Simulates a sustained flex movement to aid in building muscle tissue.

Mode 9: Tapping and Kneading – A combination of fast and slow contractions for a mixed stimulation effect.

Mode 12: Atlas – Designed for muscle confusion and athletic performance enhancement.

Indication for use

TENS:

To be used for the temporary relief of pain associated with sore and aching muscles in the shoulders, waist, back, neck, upper extremities (arms), and lower extremities (legs) caused by strain from exercise or routine household activities.

Powered Muscle Stimulator (EMS)

Designed to stimulate healthy muscles, including abdominal muscles to enhance or support overall muscle performance.

Great for:

- Stiff neck
- Sore and tight muscles
- Back pain
- Muscular performance
- Aching joints
- General pain relief

Recommended Usage Points:

Shoulders, neck, shoulder blades, along the spine, middle back, lower back, abdomen.

Arms: *biceps, triceps, forearms.*

Legs: *calves, thighs.*

Feet: *toes, heel, midfoot.*

Joints: *shoulder, elbow, wrist, hip, knee, ankle, etc.*

Unique Benefits

Quick acting – While individual results may vary, many users report noticeable relief after just 10 to 20 minutes of use.
No side effects – Since no drugs or chemicals are involved, ITouch Relief muscle stimulators are safe for long-term use. They are also safe to use alongside other medications.

Portable – Its compact, pocket-sized design allows you to enjoy treatment anywhere, at your convenience.

Cost effective – This device can save you both time and money, while enhancing your physical function and overall enjoyment of life.

Simple to use – High-tech components handle most programming automatically. The device is highly simplified for ease of use.

Be responsible – With ITouch Relief muscle stimulators, you can take charge of your own health and well-being.

Help others – Muscle tension affects nearly everyone. Now, you can help others help themselves. Whether caused by work, activity, or underlying issues—chronic or acute—muscle aches and pains are universal. Young or old, everyone needs relief and will appreciate your gift of better health.

Always available – It's ready when you need it—24 hours a day, every day. You can use it as often as needed.

Extremely effective – An advanced microchip generates multiple preset electrical, therapeutic waveforms through the electrodes. These signals deliver deep tissue pain relief and stimulate various muscle responses—helping to significantly reduce muscle tension.

Amazingly effective!

Regular Application Principles

- **Find the exact pain point:** The electrode pads should be placed directly on the area that is usually the most painful.
- **Intensity:** Gradually increase the intensity to the highest level you can comfortably tolerate.
- **Application duration:** Use the device 2–3 times per day, for 40–60 minutes each session, over the course of 10 days for optimal results.
- **There are two methods for pad placement:** twin mode and opposed mode.

Figure 1: **Twin Mode**

Fix the two pads at the top and bottom or both ends on the same side of the human body



Figure 2: **Opposed Mode**

Fix the two pads respectively on the opposite sides of the application position, as shown below

Regular Application Methods

• Pain in the Trunk

1. Probability Pain Points in the Trunk

According to Traditional Chinese Medicine (TCM), the most sensitive pressure pain point is the key target area for effective relief. Refer to Fig.1 to help identify these points. Place one electrode pad directly on the pressure pain point. Apply the second pad either symmetrically opposite the first one or close to it; but ensure the pads do not overlap. If there are multiple distinct pressure pain points, treat them one at a time.

2. Recommended Pad Positions for application on the neck, shoulders, and back (See Fig.3)

Figure 3 Use Mode 1 or 2 for 30 minutes, followed by Mode 3, 4, 5 or 6 for 30 minutes by turns for A, B, and C. For areas with significant pain, extend the treatment time and use a relatively higher intensity.


3. Recommended Pad Placement for the Waist (See Fig.4)

Figure 4 Use Mode 1 or 2 for 30 minutes, then Mode 3 or 4 for 20 minutes, or Mode 7 or 8 for 30 minutes. Light waist exercises after the session are advisable to enhance the treatment effect.

• Pain in the Joints and Limbs

Joint pain, sports injuries, and soft tissue injuries can lead to aches in areas such as the shoulder, elbow, knee, wrist, ankle, etc.

1 Position of Pads

For joint pain, apply two pads on opposite (symmetrical) sides of the affected joint. Refer to Fig.2 (opposed mode) for guidance.

2 Mode and Duration

Use Mode 1 or 2 for 30 minutes, followed by Mode 5 or 6 for 30 minutes. Apply the treatment 2–3 times daily for best results.

Pain in the muscles and tendons of the limbs is often caused by sports injuries, overexertion, or other factors. Symptoms may include muscle pain, swelling, and spasms.

1 Position of Pads


Apply one pad to each end of the affected muscle, or place one pad on the muscle and the other on the associated tendon. You may slightly adjust pad positions based on comfort and response.

2 Mode and Duration


Use Mode 3 or 4 for 20 minutes, followed by Mode 5 or 6 for 20 minutes. Consult your doctor before using the device if the injury involves a muscle or tendon tear.

Electrode Placement Variations


• Shoulders


• Hands


• Forearm


• Hip


• Feet

• Waist and Leg

• Knees

Troubleshooting

Problem	Cause	Solution
The unit is not turning on or the screen is not displaying anything.	The battery is fully depleted.	Ensure the unit is fully charged before use.
The unit does not appear to be charging.	The charging cable is not properly connected.	Confirm the USB cable is firmly connected to the unit on one end and securely plugged into the AC adapter and electrical outlet on the other. To check for a faulty AC adapter, try connecting the USB cable from the unit to a computer's USB port.
The unit turns on, but the screen is dim.	The battery is low.	Recharge the battery.
The unit is ON, but no sensation is felt from the electrode pads.	Both pads are not firmly on the skin / There is a loose connection / The intensity level needs to be increased.	Ensure you are using 2 pads per lead and they are in close proximity to each other, i.e. not on different limbs or too far apart - Make sure both pads are firmly pressed to the skin-Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.

Problem	Cause	Solution
Adhesive electrode pads do not stick to the skin, even after the surface has been cleaned and moistened.	Adhesive gel pads need to be replaced.	Replace the set of adhesive gel pads. New pads can be purchased from your local distributor or from www.itouchrelief.com .
One pad feels stronger than the other.	This is normal. Different areas will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During use, the skin feels a painful burning sensation or the stimulation becomes weaker.	Adhesive gel pads are not adhering firmly to the skin or the gel pads are too dry.	Apply a few drops of water to your fingers and rub them on both pads to help maintain the electrodes' adhesiveness. Ensure the pads are pressed firmly to the skin during use.
Did not provide sufficient relief.	Not using it for a sufficient duration or frequency.	Use your MD-9 for 20-60 minutes per session, 3-6 times a day. It is recommended to use it while lying down.
Only one port on the unit works.	Damaged port, or pads or lead wire not conducting properly.	Identify the issue by testing the problematic port with a different lead and set of pads that are known to be working.

Contraindications

- DO NOT use this device if you have any of the following medical conditions:
 1. An implanted cardiac pacemaker, defibrillator, or other implanted metallic or electronic device
 2. Undiagnosed pain syndromes
 3. Diagnosed with cancer
 4. Pregnancy
 5. Acute trauma or recent surgical procedure (within the past six months)
 6. Cardiac problems or heart disease
 7. Epilepsy
 8. Painful and/or atrophied muscles
 9. Abdominal or inguinal hernia
 10. Limited range of motion in skeletal joints
 11. Blood circulatory problems
- The MD-9 unit is intended for use by ADULTS and ONLY on healthy muscles. NEVER use this device on muscles that are atrophied, painful, experiencing spasms, or located on a limb with painful or affected joints.
- The MD-9 unit is not intended for treating any medical condition or disease, nor is it intended for physiotherapy or muscle rehabilitation. It is not to be used on any muscle that is injured or diseased.

- Do not attempt to use the MD-9 unit for muscle reeducation, prevention of muscle atrophy or spasms, improving range of motion, or treating blood flow deficiencies/venous thrombosis.

Warnings

- 1 Stimulation should not be applied over the carotid sinus nerves, especially in individuals with known sensitivity to the carotid sinus reflex.
- 2 Stimulation should not be applied transthoracically, as introducing electrical current into the heart may cause cardiac arrhythmias.
- 3 Stimulation should not be applied transcerebrally.
- 4 Apply electrodes only to normal, intact, and clean skin. Do not apply electrodes over open wounds or over swollen, infected, or inflamed areas or skin eruptions such as phlebitis, thrombophlebitis, varicose veins, etc.
- 5 DO NOT APPLY STIMULATION:
 - Over the frontal area of the neck (near the carotid sinus nerves).
 - Over the neck or mouth, as severe spasms of the laryngeal and pharyngeal muscles may occur, potentially causing difficulty in breathing. Stimulation in these areas may also adversely affect heart rhythm or blood pressure.
 - Transcerebrally.
 - Over swollen, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins, etc.).

- 23 -

- Across the chest, as the introduction of electrical current may be fatal. Consult your physician before using this device, as it may cause lethal heart rhythm disturbances in susceptible individuals.
 - Over or near cancerous lesions.
 - If you are epileptic.
 - After experiencing acute trauma or fracture.
 - Following recent surgery.
 - If you have a hernia (abdominal or inguinal).
 - To the frontal, laryngeal, and temporal regions of the neck.
- 6 Never use the MD-9 unit while driving, operating machinery, or during activities where involuntary muscle contractions could pose a risk to yourself or others.
 - 7 The effects of brain stimulation are unknown. Do not apply stimulation across the head or place electrodes on opposite sides of the head.
 - 8 Do not use the device in the bath or shower.
 - 9 Individuals with suspected heart problems or epilepsy should seek appropriate medical advice.
 - 10 Never use the MD-9 unit while sleeping.
 - 11 Never immerse the MD-9 unit in any liquid.

- 24 -

Precautions

- This MD-9 unit should be used only with the leads, electrodes, and accessories supplied by the manufacturer.
- Do not use the MD-9 unit while driving, operating machinery, or during any activity where involuntary muscle contractions may pose a risk of injury.

Adverse Reactions

- You may experience skin irritation or burns beneath the electrodes applied to your skin.
- You may experience headaches or other painful sensations during or after using electrical stimulation near the eyes, head, or face.
- Stop using the MD-9 unit and consult your physician if you experience any adverse reactions.

Location of Stimulation

Stimulation near metal: Do not use the stimulator near metal. This may increase discomfort or cause burns. Remove jewelry, piercings, belt clips, or any removable metal objects near the stimulation area.

Stimulation of the eyes: Do not apply stimulation directly to the eyes.

Stimulation of the neck or mouth: Do not use the stimulator on the neck (especially over the carotid sinus) or mouth. This may cause severe muscle spasms, difficulty breathing, and adverse effects on heart rhythm or blood pressure.

- 25 -

Stimulation of Healthy Skin: Use the stimulator only on healthy, intact skin. Stimulation may aggravate conditions involving lesions or skin diseases.

Stimulation Across the Torso: Do not apply stimulation across the torso (e.g., one electrode on the chest and one on the back). The electrical current passing through the torso can result in potentially fatal disruptions to heart rhythm.

Cranial Stimulation: The effects of stimulation on the brain are unknown. Do not apply stimulation to the head, and never place electrodes on opposite sides of the head.

Conformity to Safety Standards

EMC Statement

IEC 60601-1-2:2014

Conformity to MDD Requirements

IEC 60601-1:2005 + AMD1:2012 / ANSI/AAMI ES60601-1:2005 + A1:2012

IEC 60601-1-11:2015

IEC 60601-2-10:2012

- 26 -

Warranty Information

- This warranty covers defects in materials and workmanship under normal use for a period of 2 years from the date of retail purchase by the original end-user ("Warranty Period"). If a defect arises and a valid claim is received within this period, ITouchRelief will either repair the unit or replace it with a new product or one made from new or serviceable used parts.
- To register your product warranty, visit www.itouchrelief.com, go to the Warranties menu, and select the appropriate option (Standard 2-Year or Lifetime VIP).
- Alternatively, you may complete and mail the warranty registration card included in your package to:

Customer Service
I Touch Relief
2250 nw 114th ave,
Miami FL 33192
U.S.A.

For more information or to request repairs outside the warranty period, visit www.itouchrelief.com.

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