Fore more information, please visit our website at:
www.itouchrelief.com

True Touch

User Manual

EN

ES

User Manual

English

Contents

1	About Your I Touch Relief True Touch	2	5	Maintenance and Troubleshooting	22
1.1	Product Description	2	5.1	Maintenance, Storage & Disposal	22
1.2	Overview	3	5.2	Troubleshooting	23
1.3	How to Charge the Remote and Receivers	5	6	Technical & Legal Information	25
2	Electrode Information	6	6.1	Technical Specifications	25
2.1	Contiguous Placement	6	6.2	Conformity to Safety Standards	26
2.2	Electrode Placement	7	6.3	FCC Warning Statement	27
2.3	Application of Re-Usable Self Adhesive Electrodes	7	7	Support & Warranty	28
3	Using the Device	8	7.1	Warranty	28
3.1	Instructions for Use	9	7.2	Contact Information	29
3.2	Mode Descriptions	11			
3.3	Pulse Parameters	13			
3.4	Regular Application Principles	14			
3.5	Regular Application Methods	15			
4	Safety and Warnings	18			
4.1	Warnings	18			
4.2	Precautions	19			
4.3	Contraindications	20			
4.4	Adverse Reactions	21			
4.5	Tips for Skin Care	21			

Before operation, please read this user manual carefully and ensure you understand all instructions.

1 About Your I Touch Relief True Touch

1.1 Product Description

The I Touch Relief True Touch is a battery-powered wireless pulse generator designed to help relieve pain by sending gentle electrical impulses through the skin to stimulate nerves.

- This device works using two main parts:
- A Touch Screen Remote that controls the system.
 - A Receiver Disk with attached snap-on electrodes that stick directly to your skin.

Using the touch screen remote, you can easily adjust the intensity and choose different therapy modes. The receiver with the electrodes attaches firmly to your body and stays in place during use. The LCD display clearly shows which mode you’re using and the current intensity level, so you always know what’s happening. Refer to the figures in this manual to learn about the components and controls of your I Touch Relief system, and follow the simple steps to get the best possible pain relief.

Important

Because the I Touch Relief True Touch is completely wireless (no lead wires), only use I Touch Relief snap-on electrodes with this device.

1.2 Overview

What's Included in the Box

Your I Touch Relief True Touch kit includes:

- 1 Touch Screen Remote Control
- 2 Wireless Receivers (labeled 1 and 2)
- 1 Large Lower Back Electrode Pad (4 x 9 inches – rectangular shape, ideal for back pain relief)
- 1 Set of Extra-Large Adhesive Electrode Pads (3 $\frac{3}{4}$ x 1 $\frac{7}{8}$ inches – rectangular shape)
- 1 Set of Large Adhesive Electrode Pads (2.4 x 2.4 inches – round shape)
- 2 Electrode Wires (8 inches each, for optional wired use)
- 2 USB Charging Cables
- 1 AC Adapter
- User Manual


In case of any discrepancy between the physical product and the images shown, please refer to the actual product.


1.3 How to Charge the Remote and Receivers

Before using your I Touch Relief True Touch system for the first time, both the Remote and the Receivers may need to be charged. To maintain battery performance, charge the device every 3–6 months if not in regular use. Charging can be done using a wall outlet or a computer USB Port.

Using a Wall Outlet

Using a Computer USB Port

Important Notes

- **Remote:** Charging is complete when the battery icon shows full.
- **Receiver:** While charging, the light is purple. Charging is complete when the light turns blue.
- The AC adapter acts as a separate power source for the device. It is not part of the core medical equipment.
- The device cannot be used while charging.

2 Electrode Information

For safe and effective use of your I Touch Relief True Touch system, follow the instructions below for electrode application, placement, and care.

Important Usage Tips

- Always follow the instructions provided in the electrode pad packaging to ensure proper use and to prevent skin irritation.
- Use only I Touch Relief snap-on electrodes with this wireless system.
- Replace electrode pads when they no longer stick properly.

2.1 Contiguous Placement

This is the most common technique. It involves placing the electrodes on either side of the pain area, so the electrical current flows through or around the pain zone. For single-channel use, place one pad on each side of the painful area; especially useful for deep tissue pain on arms or legs. For example, you can place one pad on the front (anterior) and one on the back (posterior) of the affected limb to maximize stimulation.

2.2 Electrode Placement

Proper electrode placement is key to getting the best results from TENS therapy. Everyone's body responds differently, so don't hesitate to experiment with different placements to find what works best for you.

Tip: Ask for help if you have trouble reaching the area you want to treat.

2.3 Application of Re-Usable Self Adhesive Electrodes

Application Steps

1. Clean and dry the skin thoroughly with soap and water.
2. Remove the plastic film from the electrode pad.
3. Attach the pad (connected to the receiver) firmly onto the skin; directly over or near the area of pain.

Care & Storage

1. After use, store the electrodes in their resealable bag in a cool, dry place.
2. To refresh stickiness: Apply a few drops of cold water to the adhesive side and let it air dry (sticky side up).
Note: Do not oversaturate; too much water can reduce stickiness.

Safety Notes

1. Do not apply electrodes to broken, irritated, or inflamed skin.
2. Discard electrodes when they no longer adhere well.
3. Electrodes are for single-patient use only.
4. If you experience skin irritation, stop use immediately and consult a healthcare provider.

3

Using the Device

Follow these simple steps to get started with your I Touch Relief True Touch system for fast and effective pain relief.


Step-by-Step Instructions

1 Attach the Electrodes

Snap the electrode pads onto the round Receiver. Make sure both snap connectors on the Receiver are securely connected.

2 Apply to Skin

Peel off the protective film from the adhesive pads and place them directly on clean, dry skin over or near the area of pain.

3 Turn On the Devices

Remote: Press and hold the ON/OFF button until the backlight comes on.
Receiver: Press and hold the power button until the blue indicator light appears.

4 Select a Channel


On the Remote, tap the "Channel" area to choose from channel 1 - 8. Channel number will appear on the LCD screen.

5 Choose a Mode

Tap one of the 15 mode icons on the LCD screen to select your preferred therapy mode. This only works if both the Remote and Receiver are powered on.

6 Adjust the Intensity

Tap the + or – buttons on the Remote to increase or decrease the stimulation level. Adjust until the sensation feels comfortable but effective. You'll need to readjust the intensity each time you change the mode.

3.1 Instructions for Use

Remote Functions

Your Remote Control features a full touchscreen that allows you to adjust all settings quickly and easily.

- **+ / – Icons:** Increase or decrease intensity levels.
- **Channel Icon:** Choose which Receiver (1–8) you want to control.
- **15 Mode Icons:** Select from 15 pre-programmed therapy modes.
- **Default Session Time:** Each session lasts 1 hour by default. The Remote and Receiver(s) will automatically power off after 60 minutes. You can also manually turn off the device at any time.

Powering On/Off

- **Remote:** Press and hold the ON/OFF button (located on the top of the Remote) until the screen lights up. If the screen is blank, the unit is either off or needs to be recharged.
- **Receiver(s):** Press and hold the power button on the Receiver until the blue light comes on. To turn it off, press the power button again. If the light turns red, the Receiver needs to be recharged.

Lock & Unlock Function

- To conserve battery, the Remote screen will automatically turn off after 30 seconds.
- To unlock / wake up the screen, simply press the ON/OFF button once to turn the backlight on again.

Controlling Receivers & Channels

- You can control up to 6 wireless Receivers. Receiver 1 = Channel 1, Receiver 2 = Channel 2, and Up to Receiver 8 = Channel 8.
- To switch Channels:
 1. Tap the “Channel” icon on the touchscreen.
 2. Select a channel (1 to 8) that corresponds to the connected Receiver.
 3. Choose a mode by tapping one of the 15 mode icons.
 4. Adjust the intensity using the + or – buttons.
- To add more Receivers, you'll need to sync them first. See the Troubleshooting section for instructions.

Electrode Connection Tip

Each round Receiver must be connected to **1 large electrode (with 2 snap connectors)** or **2 small electrodes (each with 1 snap connector)**. Using only one small electrode will not work — the Receiver requires both connectors to be active.

Intensity Settings

- Located just below the mode icons on the Remote screen.
- Tap + to increase and – to decrease intensity.
- There are 20 intensity levels available for each channel.

Selecting a Therapy Mode

The I Touch Relief True Touch comes with 15 pre-set modes to choose from. Just tap any mode icon on the screen to instantly activate that setting.

3.2 Mode Descriptions

I Touch Relief True Touch offers 15 pre-programmed therapy modes designed for either TENS (Pain Relief) or EMS (Muscle Stimulation). Below is a breakdown of each mode and its intended purpose:

TENS Modes – For Pain Relief & Circulation Support

These modes are best for managing pain, improving blood flow, and reducing inflammation:

- **Mode 3: Tapping**
Mimics a second heartbeat. Gently taps the muscles to boost blood flow and oxygen delivery, while helping flush out inflammation and swelling from the target area.
- **Mode 4: Cupping**
Feels like a muscle being grasped, stretched, and released. Helps lengthen and loosen tight muscles, improving flexibility and mobility.
- **Mode 6: Shuffle**
A combination pattern, great for post-workout recovery. Provides a balanced mix of pulses for gentle muscle relaxation.
- **Mode 7: Ramping Wave**
Starts with the tapping rhythm of Mode 3 and gradually increases in intensity, building to a strong contraction before releasing. This mode can help boost endurance and resistance to muscle fatigue.
- **Mode 10: Reflexology**
Specifically designed for use with electrotherapy shoes. Simulates a reflexology session through socks or shoes using various rhythmic patterns to stimulate foot pressure points.
- **Mode 11: Daily Pain Relief**
Targets general aches and pains in the arms and legs caused by everyday strain or overuse.

EMS Modes – For Muscle Training & Recovery

These modes help tone, stimulate, and rebuild muscle function:

- **Mode 1: Kneading**
Simulates the feel of a real massage using thumb-like pressure and a vibrating motion. Helps relieve muscle tension and restore flexibility.
- **Mode 2: Acupressure**
A rapid, focused pulse that targets knots or tight spots in muscles. Helps break up adhesions and release built-up lactic acid.
- **Mode 5: Hot Stone Massage**
A relaxation mode that mimics a rolling manual massage, helping the soft tissues to relax and release tension.
- **Mode 7: Contract and Hold**
Provides the longest muscle contraction of all modes. Ideal for simulating flexing exercises to build muscle tone and strength.
- **Mode 9: Tapping & Kneading Combo**
Blends fast and slow contractions to create a more dynamic, engaging stimulation — useful for overall muscle engagement.
- **Mode 12: Atlas (Muscle Confusion Mode)**
Designed for athletic enhancement. Alternates stimulation patterns to promote muscle growth and adaptation through muscle confusion.

3.3 Pulse Parameters

True Touch	TABLE: PULSE parameters				
Rated Supply Voltage (V)		DC 3.7 V			
Load Resistance (Ω)		1000 Ω			
PULSE duration(μ s)	PULSE repetition frequencies (Hz)	PULSE repetition amplitudes (V)	DC component (V)	Deviation (%)	Remarks
100	53	45~85	--	± 20	Mode 1
100	6	45~85	--	± 20	Mode 2
100	1.3	45~85	--	± 20	Mode 3
100	53	45~85	--	± 20	Mode 4
100	63	45~85	--	± 20	Mode 5
100	55	45~85	--	± 20	Mode 6
100	53	45~85	--	± 20	Mode 7.1; 7.2
100	<70	45~85	--	± 20	Mode 8.1; 8.2; 8.3; 8.4;
100	<60	45~85	--	± 20	Mode 9.1; 9.2; 9.3;
100	<48	45~85	--	± 20	Mode 10.1; 10.2;
100	8	45~85	--	± 20	Mode 11
100	<60	45~85	--	± 20	Mode 12.1; 12.2; 12.3; 12.4;

3.4 Regular Application Principles

- **Find the exact pain point:** The electrode pads should be placed directly on the area that is usually the most painful.
- **Intensity:** Gradually increase the intensity to the highest level you can comfortably tolerate.
- **Application duration:** Use the device 2–3 times per day, for 40–60 minutes each session, over the course of 10 days for optimal results.
- **There are two methods for pad placement:** twin mode and opposed mode.

Figure 1: **Twin Mode**
Fix the two pads at the top and bottom or both ends on the same side of the human body



Figure 2: **Opposed Mode**
Fix the two pads respectively on the opposite sides of the application position, as shown below

3.5 Regular Application Methods

Pain in the Trunk

1. **Probability Pain Points in the Trunk:** According to Traditional Chinese Medicine (TCM), the most sensitive pressure pain point is the key target area for effective relief. Refer to Fig.1 to help identify these points. Place one electrode pad directly on the pressure pain point. Apply the second pad either symmetrically opposite the first one or close to it; but ensure the pads do not overlap. If there are multiple distinct pressure pain points, treat them one at a time.
2. **Recommended Pad Positions** for application on the neck, shoulders, and back (See Fig.3). Use Mode 1 for 30 minutes, followed by Mode 2 or Mode 3 for 30 minutes, for A, B, and C. For areas with significant pain, extend the treatment time and use a relatively higher intensity.
3. **Recommended Pad Placement for the Waist** (See Figure 4). Use Mode 1 for 30 minutes, then Mode 2 for 20 minutes, or Mode 7, 8, 9 for 30 minutes. Light waist exercises after the session are advisable to enhance the treatment effect.


Figure 3


Figure 4

Pain in the Joints and Limbs


Joint pain, sports injuries, and soft tissue injuries can lead to aches in areas such as the shoulder, elbow, knee, wrist, ankle, etc.

1. Position of Pads

For joint pain, apply two pads on opposite (symmetrical) sides of the affected joint. Refer to Fig.2 (opposed mode) for guidance.

2. Mode and Duration

Use Mode 1 for 30 minutes, followed by Mode 3 for another 30 minutes. Apply the treatment 2–3 times daily for best results.

Pain in the muscles and tendons of the limbs is often caused by sports injuries, overexertion, or other factors. Symptoms may include muscle pain, swelling, and spasms.

1. Position of Pads


Apply one pad to each end of the affected muscle, or place one pad on the muscle and the other on the associated tendon. You may slightly adjust pad positions based on comfort and response.

2. Mode and Duration


Use Mode 2 for 20 minutes, followed by Mode 3 for another 20 minutes. Consult your doctor before using the device if the injury involves a muscle or tendon tear.

Electrode Placement Variations


• Shoulders


• Hands


• Forearm


• Hip


• Feet

• Waist and Leg

• Knees

4

Safety and Warnings

4.1 Warnings

1. Do not apply stimulation over the carotid sinus nerves, especially in patients sensitive to the carotid sinus reflex.
2. Do not apply stimulation across the chest (transthoracically); electrical current near the heart may cause cardiac arrhythmias.
3. Do not apply stimulation across the brain (transcerebrally).
4. Do not apply electrodes to:
 - Open wounds
 - Swollen, infected, or inflamed skin
 - Skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins)
 - Over or near cancerous lesions
 - Over the frontal, laryngeal, or temporal regions of the neck
 - Near the mouth or throat (due to risk of spasms, difficulty breathing, or interference with heart rhythm)
 - Across the head or on opposite sides of the head
 - If you have a hernia (abdominal or inguinal)
 - If you are epileptic or have suspected heart problems
 - Following recent surgery or after acute trauma or fracture
 - While bathing or showering

Consult your physician before using this device if you have any medical condition.
5. Do not use the device while driving, operating machinery, or performing tasks where involuntary muscle movement could be dangerous.

6. Do not use the device while sleeping.
7. Never immerse the device in water or any liquid.
8. Do not use the device in the bath or shower.
9. Do not attempt to modify the device.
10. Only use accessories, pads, and parts described in this manual.
11. Do not connect the device to any equipment not specified in this manual.
12. Do not use the device simultaneously with high-frequency surgical equipment — this may cause burns or damage to the device.
13. Keep the device at least 1 meter away from shortwave or microwave therapy equipment, as it may cause interference or instability in the stimulation output.
14. Applying electrodes near the chest may increase the risk of cardiac fibrillation.
15. TENS is a pain management tool — not a cure. It does not replace pain medications or other treatments.
16. TENS works by suppressing pain signals, which are sometimes important warnings of other issues.
17. The effectiveness of TENS therapy depends on proper patient evaluation and selection by a qualified medical professional.

4.2 Precautions

1. The I Touch Relief True Touch is designed for over-the-counter (OTC) use. You do not need a prescription to operate it.
2. Use this device only with the electrode pads and accessories provided or approved by I Touch Relief. Using unapproved accessories may affect performance or safety.
3. Do not use this device while driving, operating machinery, or during any activity where involuntary muscle contractions could put you or others at risk of injury.

4.3 Contraindications

1. Do NOT use the I Touch Relief True Touch if you have any of the following medical conditions or concerns:
 - An implanted cardiac device such as a pacemaker, defibrillator, or any other metallic or electronic implant
 - Undiagnosed pain syndromes
 - A diagnosis of cancer
 - Pregnancy
 - Recent trauma or surgery (within the past 6 months)
 - Heart disease or any cardiac condition
 - Epilepsy
 - Painful or atrophied muscles
 - An abdominal or inguinal hernia
 - Limited range of motion in joints
 - Circulatory problems or vascular conditions
2. The I Touch Relief True Touch is intended for use by adults only, and only on healthy muscle tissue.
3. Never use this device on:
 - Muscles that are injured, diseased, atrophied, painful, or prone to spasms
 - Limbs with painful or impaired joints
4. This device is not intended for:
 - Medical diagnosis or treatment of any disease

- Physiotherapy or rehabilitation
- Muscle re-education
- Preventing muscle atrophy
- Managing muscle spasms
- Improving range of motion
- Stimulating blood flow or treating venous thrombosis.

4.4 Adverse Reactions

While using the I Touch Relief True Touch, some users may experience the following side effects:

- Skin irritation or burns under the electrode pads
- Headaches or discomfort during or after stimulation, especially if used near the eyes, face, or head

If you notice any adverse reactions, stop using the device immediately and consult your physician or healthcare provider.

4.5 Tips for Skin Care

To help prevent skin irritation while using the I Touch Relief True Touch 6-12, especially if you have sensitive skin, follow these simple guidelines:

1. Clean your skin using mild soap and water where electrodes will be placed. Rinse thoroughly and dry completely.
2. Trim excess hair with scissors (do not shave) to help electrodes stick better and stay comfortable.
3. Use recommended skin prep (if advised by your healthcare provider) and allow it to dry before placing electrodes.
4. Avoid stretching the skin. Apply electrodes from the center outward without pulling them tight.
5. Tape the lead wires to your skin in a loose loop to reduce tension on the electrode pads.

6. Remove electrodes gently, pulling in the direction of hair growth.
7. Use a moisturizing lotion on the area during breaks between sessions to keep skin healthy.
8. Do not apply electrodes to irritated, broken, or damaged skin.

5 Maintenance and Troubleshooting

5.1 Maintenance, Storage & Disposal

Maintenance

- Clean your skin before each use—make sure it's free from dirt, oil, or lotion.
- Optional: Spray I Touch Relief conductor spray on your fingertips and gently rub it on both electrode pads. This helps maintain the pads' stickiness.
- Turn off the device after use.
- Remove the electrodes from your skin carefully.
- Store the electrodes properly by placing the clear protective films back on them or sticking them to the electrode pad holder.
- Put everything back in the box for safekeeping.

Storage

- Store it in a cool, dry place—away from heat, moisture, and direct sunlight.
- Ensure the storage area is well-ventilated.
- Do not place heavy items on top of the device.

Disposal

- The I Touch Relief True Touch unit must be disposed of in accordance with electronic waste regulations.
- Electrode pads are disposable and should be replaced regularly; especially when they start to lose stickiness. Dispose of used pads properly.

Note: The device has an expected service life of 3 years, while the electrodes are expected to last for 2 years.

5.2 Troubleshooting

Problem	Cause	Solution
The remote is not turning on, or the screen remains blank.	The battery is fully depleted.	Make sure the REMOTE has been fully charged.
The indicator light on the RECEIVER is red	The battery is low.	Charge the RECEIVER soon to prevent it from turning off automatically.
Adhesive electrode pads do not stick to skin even after cleaning and moistening	The adhesive gel pads need to be replaced.	Replace them with new pads, available from your local distributor or at www.itouchrelief.com .
The REMOTE and/or RECEIVER are not charging	The charging cable may not be fully plugged in.	Make sure the USB cable is securely connected to the REMOTE/RECEIVER on one end and to the AC adapter plugged into an electrical outlet on the other end.

Problem	Cause	Solution
During use, the skin feels a painful burning sensation or the stimulation becomes weakened	The adhesive gel pads may not be adhering firmly to the skin or may be too dry.	Apply a few drops of water to your fingers and rub them on both pads to help restore adhesiveness. Ensure the pads are pressed firmly against the skin during application.
When pressing the Mode icons or “+”/“–” (intensity adjustors) on the REMOTE, there is no response	The RECEIVER is not turned on or not synced correctly with the REMOTE.	Make sure the REMOTE and RECEIVERS are turned on, and the channel number on the REMOTE matches the RECEIVER number. If there is still no response, re-sync them as follows: Turn on the REMOTE. Within 10 seconds, double-click the REMOTE’s ON/OFF button. Both the “Channel 1” icon and its border will flash, indicating sync mode. Touch the “Channel” area and select the desired sync channel. Double-click the RECEIVER’s power button. Once synced, the “Channel” number will automatically advance to the next, and the indicator light will stop flashing. Exit sync mode by touching the “Channel” area until it returns to “Channel 1,” or restart the REMOTE. Note: Sync must be completed within 10 seconds after turning on the REMOTE.
The REMOTE cannot control the RECEIVER properly	There is no electrical signal transmission between them.	Ensure the distance between the REMOTE and RECEIVER is less than 8 meters.

6

Technical & Legal Information

6.1

Technical Specifications

- **Channels:**
8 channels with independent intensity control
- **Power Supply:**
Remote/Receiver: DC 3.7V
- **Dimensions:**
Remote: 128.2 mm × 65.2 mm × 10.85 mm
Receiver: Ø63 mm × 14.23 mm
- **Net Weight:**
Remote: 96 g
Receiver: 25 g
- **Operating Conditions:**
Temperature: 5°C to 40°C (41°F to 104°F)
Relative humidity: 30% – 75% (non-condensing)
Atmospheric pressure: 700 – 1,060 hPa
- **Storage and Transportation Conditions:**
Temperature: -10°C to 55°C (14°F to 131°F)
Relative humidity: 10% – 90% (non-condensing)
Atmospheric pressure: 700 – 1,060 hPa

6.2

Conformity to Safety Standards

ITouch Relief declares that this device complies with the following regulatory and safety standards:

Electromagnetic Compatibility (EMC)

- IEC 60601-1-2:2014 + A1:2020

Medical Device Directive (MDD) Requirements

- IEC 60601-1:2005 + A1:2012 + A2:2020
- ANSI/AAMI ES60601-1:2005 + A1:2012 + A2:2021
- IEC 60601-1-11:2015 + A1:2020
- IEC 60601-2-10:2012 + A1:2016

FCC Compliance

- 47 CFR FCC Part 15 Subpart B
- FCC CFR Title 47 Part 15 Subpart C, Section 15.247

Wireless Coexistence

- ANSI C63.27:2017
- AAMI TIR69:2017 (R2020)
- FCC ID: 2ACD4HD-19A

6.3 FCC Warning Statement

- Any changes or modifications not expressly approved by the party responsible for compliance may void the user's authority to operate this device.
- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
 1. This device may not cause harmful interference.
 2. This device must accept any interference received, including interference that may cause undesired operation.

RF Radiation Exposure Statement

1. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter.
2. This equipment complies with RF radiation exposure limits for an uncontrolled environment.
3. The device has been evaluated to meet general RF exposure requirements and may be used in portable exposure conditions without restriction.

Part 15.105 Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception; determined by turning the equipment off and on; the user is encouraged to try one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for further assistance.

7

Support & Warranty

7.1 Warranty

All I Touch Relief TENS models include a two-year limited warranty from the date of purchase. This warranty covers the TENS/EMS stimulator REMOTE and RECEIVER only.

This warranty does not cover damage resulting from:

- Failure to follow the operating instructions
- Accidents or misuse
- Abuse or negligence
- Alteration or modification of the device
- Disassembly or repair by unauthorized personnel

To ensure your warranty remains valid, please operate the device in accordance with the provided user manual.

7.2 Contact Information

Customer Service

I Touch Relief

2250 nw 114th ave,

Miami FL 33192

U.S.A.

For more information or to request repairs outside the warranty period, visit www.itouchrelief.com.

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